Disposable Embolic Protection Device



- · Nitinol material with great hemocompatibility. Soft braided filter and round distal tip minimize the risk of trauma to the vessel.
- · Excellent axial control and kink resistance provide superior deliverability and crossability through tortuous anatomy.
- · Gold loop marker could identify the release status of the filter, which ensures the optimal positioning of the filter in the vessel.
- The unique design of guidewire casing on the filter prevents from its displacement caused by the movement of guidewire.
- · Gradient pore diameter from distal to proximal end guarantees efficient embolus capturing while maintaining continuous blood flow
- Integrated sheaths for delivery and retrieval are compatible with 5F/6F catheters and 0.014" guidewire, which could achieve rapid exchange during the procedure.

Ordering Information

Catalogue No.	Filter Diameter (mm)	Guide Wire Length (cm)	Delivery Shaft O.D.	Retrieval Shaft O.D.
TJEP03-190	3	190	3.2F	4.2F
TJEP04-190	4	190	3.2F	4.2F
TJEP05-190	5	190	3.2F	4.2F
TJEP06-190	6	190	3.2F	4.2F
TJEP07-190	7	190	3.2F	4.2F
TJEP08-190	8	190	3.2F	4.2F
TJEP03-320	3	320	3.2F	4,2F
TJEP04-320	4	320	3.2F	4.2F
TJEP05-320	5	320	3.2F	4.2F
TJEP06-320	6	320	3.2F	4.2F
TJEP07-320	7	320	3.2F	4.2F
TJEP08-320	8	320	3.2F	4.2F

Notice: Please contact the sales for detailed information before placing the order.

Embolic Coil System

C€





2D Helical Structure

3D Complex Structure

- Gradient primary diameter ranges from 0.010" to 0.014", providing different softness with better compliance.
- · Anti-unwinding structure prevents from coil unwinding which is caused by over-stretching during the embolization process.
- · Markers on proximal end of the coil and the distal end of the delivery shaft provide accurate positioning.
- Stable detachment leads to safer operation. One single coil could be detached completely within 1 second.
- This product is compatible with most of the microcatheters in the market.

Ordering Information

Minimal microcatheter inner diameter compatibility: 0.0165"

Catalogue No. (helical)	Catalogue No. (complex)	Loop Dia. (mm)	Effective Length (cm)
TJCST1,502-2D	TJCST1.502-3D	1.5	2
TJCST1.503-2D	TJCST1.503-3D	1.5	3
TJCST1,504-2D	TJCST1.504-3D	1.5	4
TJCST0201-2D	TJCST0201-3D	2	1
TJCST0202-2D	TJCST0202-3D	2	2
TJCST0203-2D	TJCST0203-3D	2	3
TJCST0204-2D	TJCST0204-3D	2	4
TJCST0206-2D	TJCST0206-3D	2	6
TJCST0208-2D	TJCST0208-3D	2	8
TJCST2.502-2D	TJCST2,502-3D	2.5	2
TJCST2.504-2D	TJCST2.504-3D	2.5	4
TJCST2,506-2D	TJCST2,506-3D	2.5	6
TJCST2.508-2D	TJCST2.508-3D	2.5	8
TJCST0304-2D	TJCST0304-3D	3	4
TJCST0306-2D	TJCST0306-3D	3	6
TJCST0308-2D	TJCST0308-3D	3	8
TJCST0310-2D	TJCST0310-3D	3	10
TJCST0312-2D	TJCST0312-3D	3	12
TJCST3.506-2D	TJCST3.506-3D	3.5	6
TJCST3.508-2D	TJCST3.508-3D	3.5	8
TJCST3,510-2D	TJCST3,510-3D	3.5	10
TJCST3.512-2D	TJCST3.512-3D	3.5	12
TJCST0404-2D	TJCST0404-3D	4	4
TJCST0406-2D	TJCST0406-3D	4	6
TJCST0408-2D	TJCST0408-3D	4	8
TJCST0410-2D	TJCST0410-3D	4	10
TJCST4.506-2D	TJCST4.506-3D	4.5	6
TJCST4.508-2D	TJCST4,508-3D	4.5	8
TJCST4.510-2D	TJCST4.510-3D	4.5	10
TJCST4,512-2D	TJCST4,512-3D	4.5	12
TJCST4.515-2D	TJCST4.515-3D	4.5	15

ı.	Catalogue No.	Catalogue No.	Loop Dia,	Effective
	(helical)	(complex)	(mm)	Length (cm)
	TJCST0509-2D	TJCST0509-3D	5	9
	TJCST0510-2D	TJCST0510-3D	5	10
	TJCST0515-2D	TJCST0515-3D	5	15
	TJCST0520-2D	TJCST0520-3D	5	20
	TJCST0610-2D	TJCST0610-3D	6	10
	TJCST0611-2D	TJCST0611-3D	6	11
	TJCST0615-2D	TJCST0615-3D	6	15
	TJCST0620-2D	TJCST0620-3D	6	20
	TJCST0715-2D	TJCST0715-3D	7	15
	TJCST0720-2D	TJCST0720-3D	7	20
	TJCST0730-2D	TJCST0730-3D	7	30
	TJCST0815-2D	TJCST0815-3D	8	15
	TJCST0820-2D	TJCST0820-3D	8	20
	TJCST0830-2D	TJCST0830-3D	8	30
	TJCST0920-2D	TJCST0920-3D	9	20
	TJCST0930-2D	TJCST0930-3D	9	30
	TJCST1020-2D	TJCST1020-3D	10	20
	TJCST1030-2D	TJCST1030-3D	10	30
	TJCST1130-2D	TJCST1130-3D	11	30
	TJCST1230-2D	TJCST1230-3D	12	30
	TJCST1330-2D	TJCST1330-3D	13	30
	TJCST1430-2D	TJCST1430-3D	14	30
	TJCST1530-2D	TJCST1530-3D	15	30
	TJCST1630-2D	TJCST1630-3D	16	30
	TJCST1830-2D	TJCST1830-3D	18	30
	TJCST2030-2D	TJCST2030-3D	20	30

Notice: Please contact the sales for detailed information before placing the order.